



Improved Bio-Agent Recovery From Simulated Environmental Surfaces with Liquid Rinse Vacuum Collection System

Detection Technologies 2008

Phoenix, AZ

Jared G. Maughan
VP of Microbiology and Operations
Microbial-Vac Systems, Inc.
Email: *Jared.Maughan@m-vac.com*



Why are We Talking About Sampling at a Detection Conference?

- All detection starts with sampling
- If you can't COLLECT it, you can't DETECT it
- If the capabilities of modern detection equipment are to be realized then we need to start with a quality sample



Calculating Efficiency: Sampling to Detection

- Sampling Efficiency (η_{Sam})

$$\eta_{\text{Sam}} = \eta_{\text{Loc}} \times \eta_{\text{Rec}}$$

- Detection Efficiency (η_{Det})

- Not all detection devices are equal
- Assume η_{Det} is 1.0 for this presentation

- Overall Process Efficiency (η_{Tot})

$$\eta_{\text{Tot}} = \eta_{\text{Sam}} \times \eta_{\text{Det}}$$



Calculating Efficiency: Sampling to Detection

- Sampling Efficiency (η_{Sam})

$$\eta_{\text{Sam}} = \eta_{\text{Loc}} \times \eta_{\text{Rec}}$$

- Detection Efficiency (η_{Det})
 - Not all detection devices are equal
 - Assume η_{Det} is 1.0 for this presentation

- Overall Process Efficiency (η_{Tot})

$$\eta_{\text{Tot}} = \eta_{\text{Sam}} \times \eta_{\text{Det}}$$

This presentation
deals with improving
this



Sampling Devices Compared

- Swab
 - Cotton tipped
 - Pre-moistened
- Sponge
 - Cellulose based
 - Pre-moistened
- Cotton Gauze
 - 2" x 2" 8 ply
 - Pre-moistened
- Microbial-Vac System (M-Vac)
 - Surface rinse
 - Vacuum collection





Overview of the M-Vac


SEC



SRS



M-Vac & Sampling Head
(MS Kit)



Vacuum

Spray



Sampling Efficiency Equation

$$\eta_{\text{Sam}} = \eta_{\text{Loc}} \times \eta_{\text{Rec}}$$




Location Efficiency (η_{Loc}) cont.

 = Contamination

 = Swab sample sites*


*Swab sample sites are 2" x 2"




4 sq ft Surface with 2" x 2" grid and 4 isolated sites of contamination



Location Efficiency (η_{Loc}) cont.




*Sponge sample sites are
4" x 4" or ~ 100 sq cm



4 sq ft Surface with 4" x 4" grid and 4 isolated sites of contamination




Location Efficiency (η_{Loc}) cont.



Method	# of samples	Total surface area covered (sq in)	η_{Loc}^*
Swab	65	260	92%
Sponge	16	256	91%
M-Vac	1	256	91%

*Assumes sampling a contaminated site would result in a positive detection result.



4 sq ft Surface with 2" x 2" grid and 4 isolated sites of contamination



Sampling Efficiency Equation

$$\eta_{\text{Sam}} = \eta_{\text{Loc}} \times \eta_{\text{Rec}}$$




Recovery Efficiency (η_{Rec})


Procedures for determining η_{Rec}

1) Spot Inoculate

- B. sub spores
- 150/surface



2) Spread inoculum over entire area





3) Allow inoculum to dry/adhere

4) Sample following standard procedures or manufactures suggestions for use

5) Elute bacteria from sampling device


- Swab: Place in 10mL, Vortex 12 times for 10 sec.
- Sponge/ Gauze: Stomach sponge and 15 mL for 2 min at 240 RPM
- M-Vac: Wrist shake 5 – 10 sec.

6) Plate on TSA plates, incubate 18 – 24 hrs, enumerate by direct plate count





Recovery of *B. subtilis* Spores from Simulated Environmental Surfaces





Reasons for Increase of η_{Rec}

Two areas of improvement using the M-Vac

- Extraction efficiency (η_{Ex})
 - Ability to remove microbes from cracks, fissures, and other surface irregularities using surface rinse solution
 - Transition from *passive* sampling to *aggressive* sampling
- Elution efficiency (η_{El})
 - M-Vac collects the bacteria in solution eliminating the need for elution
 - Overcome the sampling paradox between extraction and elution



Results of the Equation

Method	η_{Loc}	x	η_{Rec}^*	=	η_{Sam}
Swab	.02	x	.20	=	.003
Gauze	.12**	x	.23	=	.03
Sponge	.12	x	.37	=	.04
M-Vac	.91	x	.81	=	.74

*Values for η_{Rec} pertain to recovery results from stainless steel.

**Sponge value used



Conclusion

- M-Vac provides numerous advantages over historic sampling methods
 - Greater location efficiency with fewer samples to the lab (Fewer detections performed in lab)
 - Greater recovery efficiency (Improved characterization of surface contamination)
 - Increased repeatability between users
 - Increased efficiency of sample prep in the lab
 - Less exposure of lab personnel to pathogens
- Increased η_{sam} leads to increased efficiency overall ($\eta_{\text{Tot}} = \eta_{\text{Sam}} \times \eta_{\text{Det}}$)



Contributors:

Kevin Church

Kelly Black

Bruce Bradley, PhD

This material is based upon work supported by the U.S. Army

Contracting Agency – Dugway Proving Ground Directorate of Contracting under Contract No. W911S6-06-C-0004