
Comparison of the Microbial-Vac System[®] to Excision for Sampling of Beef Trimmings



Jared Maughan

VP Laboratories

Microbial-Vac Systems, Inc.

Outline

- What is the Microbial-Vac (M-Vac) and How Does it Work?
- Plant Staged Comparison of M-Vac to Excision on Trim
- Recovery of Low Level (1 – 5 CFU) *E. coli* O157:H7 from Meat Surfaces
- Summary
- Conclusions

M-Vac System

Support Equipment
Case (SEC)




M-Vac & Sampling Head
(MS Kit)



↑
Sample Collection Bottle


Surface Rinse
Solution (SRS)



How It Works



Cross section of a
piece of meat



**Spray and vacuum
action performed
simultaneously**


- Sampling head sprays sterile buffer onto surface to penetrate crevices and suspend microbes
- Solution/microbes are subsequently vacuumed up into the sample collection bottle
- Non-destructive sample acquisition
- Sample collection bottle is then removed, sealed and sent to lab

Plant-Staged Studies

- Compared the M-Vac to excision technique on beef trim at three different processing facilities
- Facilities agreed to allow quantification of APC, TC, and presence/absence of generic *E. coli*
- Conducted on the fabrication floor

Procedures for Plant-Staged Studies

- Each sample consisted of 1 site from each of 12 randomly selected pieces of meat




- Excision and M-Vac samples were in close proximity to each other on each piece of meat.
- Sites were 11.5 square cm per technique

Procedures for Plant-Staged Studies

(cont.)

- Twenty samples were taken at each facility



X 20



- Sampling technique for each site
 - Excision sites were cut to a thickness of $\frac{1}{4}$ - $\frac{1}{8}$ "
 - M-Vac sites were collected with 2 passes of the sampling head over each site

Procedures for Plant-Staged Studies

(cont.)

- Lab processing
 - Excision samples had 100 mL of diluent added and then hand massaged for 2 minutes
 - M-Vac samples were shaken for 5 seconds
- Each sample was plated onto Petrifilm™ for APC, TC, and +/- for generic *E. coli*
- Film was incubated at 37°C for 48 hours prior to direct plate count/visual inspection

Equality with Excision - Plant A





Coliform data points at zero are not shown on graph.

+/- for Generic <i>E. coli</i> *		
Method		
Sample	M-Vac	Excision
1	-	-
2	-	-
3	-	-
4	-	-
5	-	-
6	-	-
7	-	-
8	-	-
9	-	-
10	-	-
11	-	-
12	-	-
13	-	-
14	-	-
15	-	-
16	-	-
17	-	-
18	-	-
19	-	-
20	-	-
Total	0/20	0/20

*Negative (-) on Generic *E. coli* signifies that it was not present at detectable levels. Positive (+) signifies that it was.

Equality with Excision - Plant B





+/- for Generic *E. coli*

Sample	Method	
	M-Vac	Excision
1	-	-
2	-	-
3	+	-
4	-	-
5	-	-
6	-	+
7	-	-
8	-	-
9	-	-
10	-	-
11	-	-
12	-	-
13	-	-
14	-	+
15	+	+
16	+	+
17	+	+
18	+	+
19	+	+
20	+	-
Total	7/20	7/20

50% Lean (Samples 1-12)
90% Lean (Samples 13-20)

Equality with Excision - Plant C





+/- for Generic *E. coli*

Sample	Method	
	M-Vac	Excision
1	-	-
2	-	-
3	-	-
4	-	-
5	-	-
6	-	-
7	+	-
8	+	-
9	-	-
10	+	-
11	-	-
12	-	-
13	-	-
14	+	-
15	-	-
16	-	-
17	-	-
18	-	-
19	-	-
20	+	-
Total	5/20	0/20

90% Lean (Samples 1-10)
93% Lean (Samples 11-20)



Combined Plant-Staged Data Total Aerobic Bacteria Count



*P values were calculated using an un-paired t-test. No significant differences were observed between the two sampling techniques in any plant or in the combined data ($P < 0.05$).



Combined Plant-Staged Data Total Coliforms



*P values were calculated using an un-paired t-test. Significant differences were observed between the two sampling techniques in plant C ($P < 0.05$).

Coliform data points at zero are not shown on graph

Combined Plant-Staged Data Generic *E. coli*

Plant	Method	% of Samples <i>E. coli</i> +
Plant A n=20	M-Vac	0.0% (0/20)
	Excision	0.0% (0/20)
Plant B n=20	M-Vac	35.0% (7/20)
	Excision	35.0% (7/20)
Plant C n=20	M-Vac	25.0% 5/20
	Excision	0.0% (0/20)
Total n=60	M-Vac	20.0% (12/60)
	Excision	11.7% (7/60)


- No significant differences between techniques in any plant at the 5% probability level
- Plant C highlights the random nature of contamination on beef trimmings
- Given the high level of variability between the three plants, doing statistics in the combined data would be an unfair test

Low Level Contamination Testing


- Evaluate M-Vac equivalence to excision for recovery of low level (1 – 5 CFU) *E. coli* O157:H7 contamination
- Two modes of contamination
 - Pre-Rigor: Hide removal to Hot box
 - Post-Rigor: Down stream processing (cross contamination from piece to piece, knives, conveyor belts, etc.)
- For evaluation purposes pre and post refers to the condition of the meat at the time of inoculation
- Both modes were sampled post-rigor

Procedures for Low Level Contamination Testing – Lab Work

- Brisket & Flank sections collected from freshly slaughtered carcasses
 - Pieces chilled at least overnight for post-rigor work
- Outlined sites on each piece inoculated w/ 1 – 5 CFU *E. coli* O157:H7
- Inoculated sites held overnight at 2 – 4°C
- After refrigeration each site sampled with either M-Vac or Excision
- mEHEC broth added to samples and incubated at 42°C overnight
 - 5X mEHEC used for M-Vac samples to give a final concentration of 1X (25 mL 5X into 100mL sample)
 - 1:5 dilution by weight of excision sample to 1X mEHEC followed by 2 min. of hand massage
- Enriched samples then processed via standard procedures for BioControl GDS (> 98% sensitivity) with results reported as +/- for O157:H7



Low Level Contamination Data



No significant differences were observed in either treatment at the 5% probability level. A χ^2 value > 3.84 would have suggested a significant difference at the 5% probability level.

Summary

- The M-Vac technique was equivalent to excision on beef trim in the plant setting
 - Total Aerobic Bacteria
 - Total Coliform
 - Generic *E. coli*
- Low level testing shows M-Vac recovery rates that are not significantly different than excision regardless of the mode of contamination (pre vs. post-rigor)

Conclusions

- The M-Vac produces bacteria recoveries from meat surfaces that are not significantly different than excision
- The M-Vac is a suitable alternative sampling technique to excision for accomplishing N-60 sampling of trim leading to increased labor savings
- Greatly Increased Laboratory Efficiency
 - 92.5% less volume of media required*
 - No need for 2 min hand massage of each sample
 - Less incubator space required
- Non-destructive sampling
- Increased Personnel Safety (no knives, scalpels, etc.)

*Assumes 1.2 L of media for N-60 excision samples and 90 mL of 5X media for M-Vac N-60 samples.

Special Thanks . . .

- University of Tennessee
- USDA Meat Animal Research Center
- Texas A&M
- U.S. Army*
- MSI Lab Staff: Dr. Bruce Bradley, Kevin Church, Kelly Tesar, & Amanda Vos

For more information please
contact:

Glen Rose

Phone: (801) 523-3962

or

email: Glen.Rose@m-vac.com